



## Thermal Equipment Used at Our Facility

### ***ParaBath® Hand & Foot Wax Bath***



- Safe and reliable, the **ParaBath** provides deep, penetrating heat that helps relieve pain and relax stiff muscles and joints. Serves as a wonderful warm-up prior to exercise, massage and daily activities.
- Provides intensive moisture therapy that nourishes dry skin.
- Oversize inner tank measures 12" x 6" x 6", and accommodates hands, feet and elbows.
- Fully insulated unit has a worry-free temperature control that maintains paraffin at 125° to 134°F (52° to 57°C) for maximum therapeutic effect.
- Unit has a hospital-grade plug and may be run safely 24 hours a day so that it is always ready for use

### ***Hydrocollator® Moist Heat Packs***



- A hydrocollator pack is a silica gel filled pack that is soaked in hot water to provide prolonged moist heat.
- A hydrocollator pack is good in any situation that requires penetrating heat, and is an especially good choice for conditions such as hyperkyphosis that benefit from weight in addition to the heat – the larger packs are somewhat heavy when wet. For the same reason, a hydrocollator pack should not be used if the weight of the pack is excessive for the person's constitution or the tissue being treated.

### ***PneuGel® Cold Therapy Packs***



- Designed for post-op or post injury use.
- This unique ice pack support offers **PneuGel®** dual compartment technology. The bladder is made from a durable material that is inflatable and encapsulates a freezable gel. The outer pneumatic bladder compartment can be inflated to insure a comfortable and secure fit.
- Provides measured compression to reduce swelling and edema.
- A fabric layer next to the skin reduces potential for skin irritation and condensation buildup during "cold therapy".
- The wrap is made from a sturdy lightweight nylon fabric/foam cover.
- An easily accessible air valve and included larger hand pump allow for easy inflation and deflation. The support is held in place by hook and loop straps.

---

**339 Wellington Road South Unit 240, London, Ontario**  
**519-439-6111**                      [www.fpcLondon.com](http://www.fpcLondon.com)



## Thermal Equipment Used at Our Facility

### ***Cryocup® Ice Massage Devices***



- Ice Massage tool for the treatment of strains, sprains, muscle spasms, pain and inflammation of osteoarthritis and other musculo-skeletal conditions.
- Combines the cooling response of an ice pack with the effectiveness of ice massage.
- A simple-to-use ice applicator - reduces the time it takes for soothing cold to penetrate into muscles and begin relieving the aches and pains of muscle strains, sprains and tendonitis.

### ***Economy Hot/Cold Packs***



- Reusable hot/cold packs.
- Soft gel remains pliable when frozen.
- Hot/Cold Packs can be heated either in warm water on a conventional stove or microwaved in a container with water.
- Remains hot or cold for 20 minutes.

### ***Therapeutic Ultrasound for Deep Tissue Heating***



- The [Rich-Mar](#) Ultrasound offers fantastic flexibility in a lightweight ultrasound that will work equally well in the clinic or on the road.
- Rich-Mar's patented "Therapy Hammer" transducer that provides single transducer with both a 2cm<sup>2</sup> and a 5cm<sup>2</sup> sound head.
- Both the 2cm<sup>2</sup> and a 5cm<sup>2</sup> can output both 1 (for deep penetration) & 3 (for superficial penetration) MHz frequencies.
- Continuous (heating effect) and pulsed (healing effect) outputs.
- Adjust time, frequency, transducer size, intensity any time during the treatment.

---

**339 Wellington Road South Unit 240, London, Ontario**  
**519-439-6111**                      [www.fpcLondon.com](http://www.fpcLondon.com)